






Inhalational Anaesthesia Equipment

EQUIPMENT	TICK
Vaporizers × 2	
Sevoflurane gas* × 2	
Scavenging tubing × 2	
Easy-fil bottle adapter	
Scavenging tubing set up instructions (Available on Right Decisions – Guidelines > Resus > Aide Memoire)	<p>Scavenging tubing set up</p> <ul style="list-style-type: none"> • Attach one end of the scavenging tubing to the expiratory valve of the C circuit • Attach the other end to the scavenging system on the anaesthetic machine 

Date box made-up:

Completed by (sign & print):

Shortest expiry date:

*Available from CEPOD OPD